



# PROSTATE CANCER BIOREPOSITORY NETWORK

SOP No: 006  
RNA Extraction from Frozen Tissues

---

STANDARD OPERATING PROCEDURE	<b>SOP No. 006</b> <b>RNA Extraction from Frozen Tissue</b> <b>UW Network Site</b>
<b><i>UNAUTHORIZED COPYING PROHIBITED</i></b>	<b>Version Number: 1.0 (November 2016)</b> Number of Pages: 4
Electronic Filename: PCBN.SOP06.v1.0 RNA Extraction Frozen Tissues	



# PROSTATE CANCER BIOREPOSITORY NETWORK

SOP No: 006

RNA Extraction from Frozen Tissues

## 1. PURPOSE

To describe the procedure for the extraction of RNA from frozen tissue. NOTE: This SOP does not cover detailed procedures for handling Human Biological Materials or hazardous chemicals and it is recommended that personnel following this SOP refer to institutional safety guidelines.

**The following SOP for RNA extraction from frozen tissue is followed by the UW Network Site of the PCBN.**

## 2. RESPONSIBILITIES

**Authorized personnel** extracting RNA from frozen prostate tissue must ensure that:

- all procedures are followed correctly
- all samples are adequately coded during processing
- all documentation is completed, and accurate records maintained on all samples

## 3. HEALTH AND SAFETY

Personnel carrying out this procedure must maintain safe working practices and observe all relevant Environmental Health and Safety (EH&S) guidelines. This includes the appropriate use of Personal Protective Equipment (PPE), and the procedures for waste disposal, disinfection & spill clean-up and biosafety.

## 4. EQUIPMENT AND MATERIALS

Equipment	Materials
Centrifuge Cryostat Fume hood gentleMACS Dissociator Refrigerated Microcentrifuge Vortex	Chloroform Cryovials 75% ethanol isopropyl alcohol gentleMACS M tube Microcentrifuge tubes Phase Lock Gel Heavy 2ml tube (5 Prime) RNA STAT-60 RNase-free water Scalpel blades



## 5. PROCEDURES

### 5.1 Sample Preparation

- Keep previously frozen tissue tubes on dry ice and frozen until in RNA-Stat 60

*If starting with frozen bit of tissue:*

- Add 1.1 ml RNA STAT-60 to gentleMACS M tube
- If tissue is less than 100 mg, place directly into M tube, if tissue is more than 100 mg, remove from tube and cut piece to less than 100 mg using a new scalpel blade before placing in M tube.
- After replacing lid, ensure that all of the sample material reaches the area of the roto/stator.
- Run the gentleMACS Program RNA\_02 (this is for frozen tissue).
- Centrifuge the M Tube at 2000 RCF for 1 minute at room temperature to collect lysate at the bottom of the tube.
- Remove homogenized sample from the tube and place into a 1.5 ml tube and let incubate at room temperature for 5 minutes.

*If starting with 1.5 ml tube of frozen tissue block sections (these should be  $\leq 20$  um thick and weigh less than 100 mg):*

- Add 1.1 ml of RNA STAT-60 to tube of tissue sections.
- Vortex for 1 minute to lyse tissue sections.
- Incubate for 5 minutes at room temperature.

### 5.2 RNA Extraction

- Pre-spin Phase Lock Gel Heavy 2ml tube (5 Prime) at 12,000 RCF for 1 minute at 4°C and set aside
- Add 100 ul of nuclease-free water to above pre-spun tube
- Centrifuge samples at 12,000 RRCF for 3 minutes at 4°C to pellet cellular debris and allow possible lipid content to layer on top of sample. If lipid content is observed at the surface of the STAT-60 samples, then pipet to remove.
- Remove 1 ml of STAT-60 sample supernatant by pipetting and place in the previously pre-spun Phase Lock Gel tube. Any remaining RNA STAT-60 can be discarded.
- Add 200 ul of chloroform.
- Cap sample tubes securely. Shake tubes vigorously by hand for 15 secs and incubate at room temperature for 2-3 min.
- Centrifuge the samples at no more than 12,000 RCF for 15 min at 4°C.
- Obtain new 1.5 ml tube, label and add 600ul of isopropanol.
- Transfer the upper aqueous phase to the new 1.5 ml microcentrifuge tube containing isopropanol (this can be easily completed by pouring as the Phase Lock Gel forms a tight seal between the phases).
- Vortex 15 seconds. Allow RNA to precipitate 5 to 7 minutes at room temperature.
- Centrifuge at 12,000 RCF for 10 minutes at 4°C to pellet RNA.



# PROSTATE CANCER BIOREPOSITORY NETWORK

SOP No: 006

RNA Extraction from Frozen Tissues

---

- Remove supernatant and discard.
- Add 1 ml of cold 75% ethanol, invert tube several times to wash RNA pellet.
- Centrifuge at 7,500 RCF for 5 minutes at 4°C.
- Remove supernatant and air dry pellet at room temperature for 5-10 minutes.
- Solubilize RNA in an appropriate amount of cold nuclease-free water.