



# PROSTATE CANCER BIOREPOSITORY NETWORK

SOP No: 007  
QC of RNA

---

STANDARD OPERATING PROCEDURE	<b>SOP No. 007</b> <b>QC of RNA</b>
<b><i>UNAUTHORIZED COPYING PROHIBITED</i></b>	<b>Version Number: 2.0 (November 2016)</b> Number of Pages: 3
Electronic Filename: PCBN.SOP07.v2.0 QC of RNA	



# PROSTATE CANCER BIOREPOSITORY NETWORK

SOP No: 007  
QC of RNA

## 1. PURPOSE

To describe the procedure for the QC of RNA extracted from frozen tissue by real-time PCR. NOTE: This SOP does not cover detailed procedures for handling Human Biological Materials or hazardous chemicals and it is recommended that personnel following this SOP refer to institutional safety guidelines.

**The following SOP for QC of RNA is followed by the JHU, NYU, MSKCC and UW Network Sites of the PCBN.**

## 2. RESPONSIBILITIES

**Authorized personnel** performing QC of RNA extracted from frozen prostate tissue must ensure that:

- all procedures are followed correctly
- all samples are adequately coded during processing
- all documentation is completed, and accurate records maintained on all samples

## 3. HEALTH AND SAFETY

Personnel carrying out this procedure must maintain safe working practices and observe all relevant Environmental Health and Safety (EH&S) guidelines. This includes the appropriate use of Personal Protective Equipment (PPE), and the procedures for waste disposal, disinfection & spill clean-up and biosafety.

## 4. EQUIPMENT AND MATERIALS

Equipment	Materials
Agilent Bioanalyzer iCycler real time Thermocycler Nanodrop	96-well plates DNase I (Ambion, Catalog # AM2222) Primers as described SuperScript First Strand Synthesis System for RT-PCR (Invitrogen, Cat # 11904-018) SYBR Green PCR Master Mix (Invitrogen Catalog # 4309155)



## 5. PROCEDURES

- Quantitate RNA samples extracted using **PCBN SOP006 RNA Extraction from Frozen Tissue** by Nanodrop. Refer to the Nanodrop user guide for instructions of use.
- Obtain RIN number for samples extracted using **PCBN SOP006 RNA Extraction from Frozen Tissue** by Agilent Bioanalyzer.

### Procedures for JHU and NYU Network Sites Only:

- DNase I Treatment (Ambion, Cat # AM2222)
  - Treat 50 µl of RNA (<200 µg/ml) with 1 µL DNase I.  
Refer to standard protocol for 'Removal of Contaminating Genomic DNA from RNA' (37°C for 30 min).
  - Heat inactivate at 75°C for 10 min following addition of EDTA to a final concentration 5 mM.
- cDNA synthesized using SuperScript First Strand Synthesis System for RT-PCR (Invitrogen, Cat # 11904-018) following standard protocol for 'First-Strand Synthesis Using Random Primers'.
- QC RNA samples extracted using **PCBN SOP006 RNA Extraction from Frozen Tissue** by real-time PCR.
  - Rxn mix (20 µl): µl/rxn

cDNA	2
SYBR Green	10
F-primer	1
R-primer	1
dH2O	6
  - GAPDH primers (81 bp product)  
  
Forward: 5'-CGCTCTCTGCTCCTCCTGTT-3'  
Reverse: 5'-CCATGGTGTCTGAGCGATGT-3'
  - 18S primers (110 bp product)  
  
Forward: 5'-GATGGTAGTCGCCGTGCC-3'  
Reverse: 5'-GCCTGCTGCCTTCCTTGG-3'
  - Cycling parameters  
  
Cycle 1: ( 1X)  
Step 1: 95.0°C for 03:00  
  
Cycle 2: ( 35X)  
Step 1: 95.0°C for 00:30  
Step 2: 60.0°C for 00:30  
Step 3: 72.0°C for 00:45